

# BIOREACTOR ENGINEERING Chapter 3 Culture Kinetic Study of Continuous Fermentation

by Chew Few Ne Faculty of Chemical & Natural Resources Engineering cfne@ump.edu.my



Culture Kinetic Study of Continuous Fermentation by Chew Few Ne

## **Chapter Description**

- Topic Outcome
  - Perform calculation regarding culture kinetics of continuous fermentation
- References
  - Doran, P.M. (2013) Bioprocess Engineering Principles. Elsevier.
  - Liu, S. (2013) Bioprocess Engineering: Kinetics,
     Biosystem, Sustainability and Reactor Design. Elsevier.
  - Rao, D.G. (2010) Introduction to Biochemical Engineering. McGraw Hill.



## **Topic Outline**

- Kinetics of Cell Growth
- Kinetics of Substrate Consumption
- Kinetics of Product Formation



• Material balance around the chemostat on the cell:

Cell accumulation = cell in – cell out + cell growth – cell death

#### Please write the formula = ?

F = volumetric flow rate of medium (feed & effluent) (L/h)

 $S_0$  = substrate concentration in the feed (g/L)

$$X_0$$
 = cell centration in the feed (g/L)

$$P_0$$
 = product concentration in the feed (g/L)

S = substrate concentration in the fermenter (g/L)

X = cell concentration in the fermenter (g/L)

- P = product concentration in the fermenter (g/L)
- V = liquid volume of fermenter (L)
- $\mu$  = specific growth rate (h<sup>-1</sup>)
- $k_d$  = specific death rate (h<sup>-1</sup>)



- For single-stage chemostat, medium supply is usually sterile
- During exponential phase:  $k_d \ll \mu \rightarrow k_d \approx 0$

• So, 
$$\frac{dX}{dt} = -\frac{F}{V}X + \mu X$$

- Let dilution rate, D:  $D = \frac{\text{medium flowrate}}{\text{culture volume}} = \frac{F}{V}$
- So, Please write the formula = ?
- During steady state,  $\frac{dX}{dt} = 0$   $\therefore \mu = D$
- The equation shows that by varying the medium supply rate, the growth rate can be varied

SO Culture Kinetic Study of Continuous Fermentation by Chew Few Ne

• From Monod equation:

$$\mu = \mu_{max} \frac{S}{K_s + S}$$

• Thus,  $D = \mu_{max} \frac{S}{K_S + S}$  OR  $S = \frac{K_S D}{\mu_{max} - D}$ 



Culture Kinetic Study of Continuous Fermentation by Chew Few Ne

• Exercise 1



Culture Kinetic Study of Continuous Fermentation by Chew Few Ne

• Material balance around the chemostat on the limiting substrate:

substrate accumulation = substrate in – substrate out – substrate for growth– substrate for product formation – substrate for maintenance

Please write the formula = ?

- F = volumetric flow rate of medium (feed & effluent) (L/h)
- V = liquid volume of fermenter (L)
- $S_0$  = limiting substrate concentration in the feed (g/L)
- S = limiting substrate concentration in the fermenter (g/L)

 $Y_{X/S}$  = the theoretical yield of cell from substrate (g cell dry weight formed per g substrate consumed)

 $Y_{P/S}$  = the theoretical yield of product from substrate (g product formed per g substrate consumed)  $\mu$  = specific growth rate (h<sup>-1</sup>)

 $q_p$  = specific product formation rate (h<sup>-1</sup>).

m = the maintenance coefficient (g substrate consumed per g cell dry weight per hour)



- In general,  $mX \ll \frac{\mu X}{Y_{X/S}}$  and can be neglected
- If no product is formed, so: Please write the formula = ?
- During steady state,  $\frac{dS}{dt} = 0$  and  $\frac{F}{V} = D$

• So: 
$$X = Y_{X/S}(S_0 - S)$$

• With 
$$S = \frac{K_S D}{\mu_{max} - D} \rightarrow X = Y_{X/S} \left( S_0 - \frac{K_S D}{\mu_{max} - D} \right)$$

Culture Kinetic Study of Continuous Fermentation by Chew Few Ne

- From the graph of cell and substrate concentration as a function of dilution rate,
  - At low D, all the substrate is consumed at steady state, S ≈ 0 and X is high:
     X =  $Y_{X/S}S_0$
  - As D increases, S increases and X decreases.
- For X falls to zero rapidly,
  - This is called a washout condition.
  - The rate at which cells are removed in the outlet stream is greater than the rate of generation by growth.
  - Although  $\mu$  can be manipulated by changing D, the cell growth rate is ultimately limited to  $\mu_{\text{max}}$
  - During washout,  $D = D_c$  (critical dilution rate)



• To calculate the value of D<sub>c</sub>, from previous:

$$D = \mu_{max} \frac{S}{K_s + S}$$

• If  $D = D_c$ ,  $S = S_0$  when X = 0, and  $K_s << S_0$ 

• So,  $D_c = \mu_{max}$ 

- We must avoid washout in practical operation where  $D < D_c$
- When the D  $\approx$  D<sub>c</sub>, small changes in D cause large changes in X and S (due to fluctuation in the flow rate).



#### **Kinetics of Product Formation**

• The cell productivity of a chemostat:  $P_X = DX$ 

 $P_X$  = volumetric rate of cell production (g/L.h) D = dilution rate (h<sup>-1</sup>) X = cell concentration (g/L)

• From previous:

$$X = Y_{X/S} \left( S_0 - \frac{K_S D}{\mu_{max} - D} \right)$$

• So, 
$$P_X = DY_{X/S} \left( S_0 - \frac{K_S D}{\mu_{max} - D} \right)$$



#### **Kinetics of Product Formation**

• From the graph of cell productivity as a function of dilution rate, the dilution rate giving rise to maximum productivity,  $D_{opt}$  can be calculated by solving:  $\frac{dP_X}{dP_X} = 0$ 

• So, 
$$D_{opt} = \mu_{max} \left( 1 - \sqrt{\frac{K_S}{K_S + S_0}} \right)$$

• If  $K_s << S_0$ ,  $D_{opt} \approx \mu_{max} \approx D_c$  (Washout may occur)



#### **Kinetics of Product Formation**

• Exercise 2



Culture Kinetic Study of Continuous Fermentation by Chew Few Ne



# CREDITS

Special thanks to
Prof. Dr. Tey Beng Ti
Prof. Madya Dr. Rosfarizan binti Mohamad
Dr. Farhan binti Mohd Said

Thanks to Dr. Rozaimi bin Abu Samah for proofreading the learning contents.



Culture Kinetic Study of Continuous Fermentation by Chew Few Ne